


pCLIP_f-H2B Control Plasmid



N9218S 002141117111

N9218S

20 µg Lot: **0021411**
Store at: -20°C Exp: **11/17**



Live CHO-K1 cells transiently transfected with pCLIP_f-H2B. Cells were labeled with CLIP-Cell™ 505 (green) for 30 minutes and counterstained with Hoechst 33342 (blue).

Introduction

This control plasmid contains the gene encoding the Histone H2B protein cloned upstream of the CLIP_f coding sequence in pCLIP_f as a fusion to the N-terminus of the CLIP-tag. Histone H2B is a member of the core histones that package DNA in the nucleus. The H2B-CLIP_f fusion protein gives nuclear fluorescence when labeled with CLIP-Cell™ substrates. The full sequence for pCLIP_f-H2B can be downloaded at www.neb.com.

The CLIP-tag is a novel tool for protein research, allowing the specific, covalent attachment of virtually any molecule to a protein of interest.

The CLIP-tag is a small protein based on human O⁶-alkylguanine-DNA-alkyltransferase (AGT). CLIP-tag substrates are derivatives of benzylcytosine (BC). In the labeling reaction, the substituted benzyl group of the substrate is covalently attached to the reactive cysteine of CLIP-tag forming a stable thioether link. Although CLIP-tag is based on the same protein as SNAP-tag®, the benzylcytosine substrates form a separate class of substrates, different from the benzylguanine substrates recognized by SNAP-tag. CLIP-tag and SNAP-tag can be used for orthogonal simultaneous labeling.

pCLIP_f contains an improved version of CLIP-tag, termed CLIP_f. CLIP_f displays faster kinetics in *in vitro* labeling and fast, specific and efficient labeling in live and fixed cell applications, thereby rendering it a desired research tool for analysis of protein dynamics.



1-800-632-7799
info@neb.com
www.neb.com

There are two steps to using this system: sub-cloning and expression of the protein of interest as a CLIP_f fusion, and labeling of the fusion with the CLIP-tag substrate of choice. Expression of the CLIP_f-H2B fusion protein is described in this document. The labeling of fusion proteins with CLIP-tag substrates is described in the instructions supplied with CLIP-tag substrates.

Materials Required but not Supplied:

Tissue culture media and reagents
Mammalian cell lines
Transfection reagents
CLIP-tag substrates

Storage

pCLIP_f-H2B is supplied in TE buffer (10 mM Tris-HCl, pH 8.0, 1 mM EDTA) at a concentration of 0.5 µg/µl. Plasmid solutions can be stored at 4°C for up to one week. For long-term storage -20°C is recommended.


Expression of CLIP_f Fusions

Transient Expression

Expression of the fusion protein cloned in pCLIP_f-H2B can be achieved by transiently transfecting cells in culture with standard transfection protocols. The appropriate reagent and time to permit adequate expression must be empirically determined. pCLIP_f-H2B has performed well in stable and transient transfection of CHO-K1, COS-7, U-2 OS and NIH 3T3 cells. Note that the intensity of the fluorescence may vary depending on cell line and labeling substrate used.

Stable Expression

pCLIP_f-H2B can be transfected as described above for transient transfection or by other standard transfection methods. Twenty four to 48 hours after transfection, begin selecting mammalian cultures in 600–1,200 µg/ml G418 (geneticin) depending on the cell line. It is recommended that a kill curve be established for each cell line to determine optimal selection conditions. After 8–12 days of continuous selection, stable colonies will become visible. It is possible to use pools of stable cell populations for initial cell labeling to test for the presence of CLIP-tag expression. In addition, clonal cell lines can be isolated and characterized if desired.



Troubleshooting

Expression

In general, we have not experienced problems expressing H2B-CLIP_f from the pCLIP_f-H2B plasmid. Labeling of transfected cells with a fluorescent CLIP-Cell substrate should show strong nuclear fluorescence. In most instances, difficulties in expression can be resolved by altering the transfection protocol.



Notice to Buyer/User: The Buyer/User has a non-exclusive license to use this system or any component thereof for **RESEARCH AND DEVELOPMENT PURPOSES ONLY**. Commercial use of this system or any components thereof requires a license from New England Biolabs, Inc., 240 County Road Ipswich, MA 01938. For detailed information, see: www.neb.com/cia/legal.

The products and/or their use may be covered by one or more of the following patents and patent applications: U.S. Patent No. 7,939,284 (Methods for Using O⁶-Alkylguanine-DNA-Alkyltransferases); U.S. Patent No. 7,888,090 (Mutants of O⁶-Alkylguanine-DNA-Alkyltransferases); U.S. Patent No. 8,163,479 (Specific Substrates for O⁶-Alkylguanine-DNA-Alkyltransferases); U.S. Patent No. 8,178,314 (Pyrimidines reacting with O⁶-Alkylguanine-DNA-Alkyltransferases); PCT/EP2007/057597 (Labeling of Fusion Proteins with Synthetic Probes); EP07117800 (Drug Delivery); EP07117802 (Drug Delivery); EP07120288 (GTPase-Transient Protein Protein Interactions) These patents and patent applications are owned by Covaly, or owned by the Ecole Polytechnique Fédérale de Lausanne (EPFL) and exclusively licensed to Covaly and NEB.

For Research Use Only: No other use is authorized, including without limitation Commercial Use. Commercial Use means any and all uses of the product and derivatives by a party for monetary or other consideration including, but not limited to: (1) product manufacture; (2) providing of a service, information or data; (3) resale of the product, derivatives or components; (4) use of the product, derivatives or components by therapeutic, diagnostic or prophylactic clinical or commercial purposes. Information regarding any other license may be obtained by contacting Business Development and Licensing, New England Biolabs, Inc., 240 County Road, Ipswich, MA 01938. Tel: 978-927-5054, Fax: 978-921-1350.

FuGENE® is a registered trademark of Roche.